

Characterization of the IMRT and SBRT performance of a novel biology-guided radiotherapy (BgRT) machine using ArcCHECK

D. Zaks, M Narayanan, R. Bassalow, O. Volotskova, Y. Voronenko, D. Pal, D. Rigie, J. Burns, A. Purwar, P. Olcott, G. Kuduvali
Reflexion Medical, Hayward, California, United States

INTRODUCTION

The Reflexion™ X1 biology-guided radiotherapy (BgRT) machine consists of a rotating gantry that includes a 6MV linear accelerator, kilovoltage fan-beam CT, two 90° PET detector arcs and a megavoltage CT detector, continuously spinning at 60 rpm. It can achieve a nominal dose rate of 850 cGy/min and supports two clinical field sizes: 40x1 and 40x2 cm².

This work describes performance measurements characterizing the Reflexion X1 dose delivery accuracy of IMRT/SBRT treatments for static targets using the ArcCHECK device.

AIM


To characterize the IMRT and SBRT performance of the Reflexion X1 BgRT machine using the ArcCHECK dosimetry device

RESULTS

Table 1: Results of Static IMRT/SBRT tests using the 3%/3mm gamma criterion for dose accuracy

Target Shape	Surrounding Material	Jaw Size	Prescription Dose	Gamma pass rate
22mm sphere	Homogenous	2 cm	1000 cGy	99.4%
22mm sphere	Heterogenous	2 cm	1000 cGy	96.9%
C-shape	Homogenous	1 cm	1000 cGy	91.0%
C-shape	Homogenous	2 cm	1000 cGy	100%
22mm sphere + C-shape (2-target plan)	Homogenous	2 cm	1000 cGy	98.0%
TG119 (Prostate)§	Homogenous	2 cm	1000 cGy	97.7%
TG119 (C-shape) §	Homogenous	1 cm	200 cGy	93.4%
TG-119 (Head-Neck) §	Homogenous	2 cm	200 cGy	97.3%

§Targets and OARs based on TG-119 were simulated rather than physically constructed



METHOD




Figure 1: Images of the different custom inserts used with the ArcCHECK: a) Homogenous insert with spherical target, b) Heterogenous insert with Styrofoam mesh and spherical target and c) Homogenous insert with two targets (C-shape + spherical)

Custom Inserts for ArcCHECK

- Customized Inserts were designed to fit inside the ArcCHECK cavity
- Inserts simulated homogenous and heterogenous material, surrounding the targets.
- Inserts accommodated different target shapes: sphere and C-shape)
- Inserts that can be configured to hold up to 2 targets were used to test delivery accuracy of IMRT/SBRT plans

Dose Delivery Accuracy

- Dose delivery accuracy was evaluated using 3%/3mm Gamma criterion [1]




Figure 2: Front-side view of the ArcCHECK without an insert (top) and front-side view of the ArcCHECK with custom insert (bottom)

Static Targets

- All eight delivered IMRT/SBRT plans met the gamma index PASS criteria of at least 90% of the measurement points achieving the 3%/3 mm closeness criteria (see Table 1)

CONCLUSIONS

The performance measurements indicate that the Reflexion X1 System has the required dose delivery accuracy to deliver IMRT/SBRT treatments in a variety of test conditions.

ACKNOWLEDGEMENTS

We would like to thank Sun Nuclear Corporation for their continued support and engagement in developing a robust QA program.

REFERENCES

- [1] D. A. Low et al, "A technique for the quantitative evaluation of dose distributions", Med Phys 1998; 25: 656-661
- [2] G. A. Ezzell et al, "IMRT commissioning: Multiple institution planning and dosimetry comparisons, a report from AAPM Task Group 119", Med Phys 36; 2009: 5359-5373

CONTACT INFORMATION

Daniel Zaks - dzaks@reflexion.com

