



Imaging the effective point of measurement of a dosimeter inside a phantom

Swapna Lilly Cyriac¹, Jian Liu², Emel Calugaru², Jenghwa Chang^{2,3}

¹Department of Radiation Oncology, KIMS Cancer Care and Research Center Pvt Ltd, Kerala, India- 695029

²Department of Radiation Medicine, Northwell Health - Donald and Barbara Zucker School of Medicine at Hofstra/Northwell, Lake Success, New York- 11042, USA

³Department of Physics and Astronomy, Hofstra University, Hempstead, New York, USA



INTRODUCTION

The Leksell Gamma Knife (LGK) solid water (SW) phantom is designed for the PTW 31010 ionization chamber for reference dosimetry (Fig 1). When this phantom is used, **effective point of measurement (EPM)** of other detectors does not coincide with the EPM of PTW 31010. In this study, we developed a simple and effective method to determine the EPM of a detector inside the existing LGK-SW phantom.

AIM

To develop an **imaging method** for identifying the **EPM** of a detector inside a **LGK-SW dosimetry phantom**.

METHOD

- When the PTW-31010 is inserted into LGK-SW phantom, the EPM is at $(x,y,z) = (100,100,100)$ of Leksell Coordinate System (LCS), which might not be same for other detectors, particularly in the z-direction.
- A strategy was developed for identifying the **z-LCS coordinate** of EPM for four diode detectors (Table 1) using the **CBCT of Gamma Knife Icon unit**.
- CBCT image was acquired for each detector and the flat face/tip of the detector was identified. The z-LCS coordinate of **EPM** was then determined based on the **published distance between EPM and the flat face/tip**.
- For validation, a plan was generated to measure the output of 4-mm collimator for PTW-TN60017 at various z-locations.

RESULTS

- Fig 2 illustrates how z_{eff} was determined for the CBCT scan of PTW-60019. In Fig2, z-LCS coordinate of the detector tip was identified to be $z = 96.5$. Since EPM of this detector is 1.3 mm from the tip, the exact EPM location was $z_{\text{eff}} = 96.5 + 1.3 = 97.8$, or 2.5 mm from the EPM of PTW-31010.
- Applying the similar calculation, z-LCS coordinate was at $z_{\text{eff}} = 98.5, 98.9$, and 97.5 for respectively PTW-TN60008, TN60016, and TN60019.
- Fig 3 shows the relative output profile along the z direction, measured using the PTW-TN60017 diode. The 50% relative output (green box) occurs at $z=95.2$ mm and 100.2 mm, so the center of the measured profile was located at $z = 97.7$, only **0.1 mm difference** from $z_{\text{eff}} = 97.8$ for this diode identified using the CBCT method.

Detector	Sensitive Volume (mm ³)	Reference point from Tip (mm)	z- LCS of EPM (mm)
PTW TN60008	0.03	2.0	98.5
PTW TN60016	0.03	2.4	98.9
PTW TN60017	0.03	1.3	97.8
PTW 60019	0.004	1.0	97.5

Table 1: Characteristics of the four diode detectors used in this study. The last column lists the z - LCS of EPM when it is inserted into the SW phantom, derived using the CBCT method.




Fig 1. Setup for the verification of EPM of different diode detectors in a LGK-SW.




Fig 2. CBCT images for PTW 60019: (a) axial, (b) sagittal and (c) coronal view. As shown in (b) and (c), $z = 96.5$ (mm) is the z-LCS of the tip of detector in the superior-inferior (S-I) direction. Since the EPM is 1 mm inferior from the tip for this diode, $z_{\text{eff}} = 96.5 + 1.3 = 97.8$ (mm) is the Z LCS of EPM.




Fig 3. The 50% relative output (green box) occurs at $z = 95.2$ and 100.2 , so the center of the profile is located at $(95.2 + 100.2) / 2 = 97.7$, which is 0.1 mm different from the identified location ($z_{\text{eff}} = 97.8$, red line) of the sensitive volume for this diode.

CONCLUSIONS

- A **simple** and **effective method** was developed to determine the EPM of diode detectors inside the holder of LGK -SW phantom designed for PTW 31010 ion chamber.
- With the acquired information, **each diode can be positioned at the focal point of the machine** for accurate relative dosimetry measurements, without any modification or customization of the LGK SW phantom.

ACKNOWLEDGEMENTS

We are grateful for the financial support from the Union for International Cancer Control (UICC), Giuseppe Motta 1202 Geneva - Switzerland. We would like to thank Mr. Sherin Joseph, Dept. of Radiation Medicine, CFAM, Northwell Health, NY for all the help and support.

REFERENCES

- An International Code of Practice for Reference and Relative Dose Determination. Dosimetry of Small Static Fields Used In External Beam Radiotherapy. **Technical Reports Series No. 483**; IAEA: 2017.
- Ma L, Kjall P, Novotny J et al. A simple and effective method for validation and measurement of collimator output factors for Leksell Gamma Knife Perfexion. **Phys. Med. 2009; Biol.54:3897-3907**
- Hrsak Hrvoje, Marija Majerb, et al. Correction of measured Gamma-Knife output factors for angular dependence of diode detectors and PinPoint ionization chamber. **Physica Medica. 30(8):914-919**.

CONTACT INFORMATION

Swapna Lilly Cyriac
Department of Radiation Oncology
Kerala Institute of Medical Science (KIMS)
Trivandrum, Kerala, India- 695029.
Email : swapnalillycyriac@gmail.com

