# The importance of Monte Carlo (Geant4) physics settings for clinical proton therapy treatment plans.

Carla Winterhalter<sup>1,2</sup>, Michael Taylor<sup>1,2</sup>, David Boersma<sup>3,4</sup>, Alessio Elia<sup>4</sup>, Susanna Guatelli<sup>5</sup>, Ranald Mackay<sup>1,2</sup>, Karen Kirkby<sup>1,2</sup>, Lydia Maigne<sup>6</sup>, Vladimir Ivanchenko<sup>7,8</sup>, Andreas Resch<sup>9</sup>, David Sarrut<sup>10</sup>, Peter Sitch<sup>2</sup>, Marie Vidal<sup>11</sup>, Loïc Grevillot<sup>4</sup> & Adam Aitkenhead<sup>1,2</sup>

<sup>1</sup> Division of Cancer Sciences, University of Manchester, M13 9PL, Manchester, UK <sup>2</sup> The Christie NHS Foundation Trust, M20 4BX, Manchester, UK <sup>3</sup> ACMIT Gmbh, Viktor-Kaplan-Straße 2, A-2700 Wiener Neustadt, Austria <sup>4</sup> EBG MedAustron GmbH, Marie Curie-Straße 5, A-2700, Wiener Neustadt, Austria <sup>5</sup> Centre For Medical Radiation Physics, University of Wollongong, Australia <sup>6</sup> Laboratoire de Physique de Clermont, Campus Universitaire des Cézeaux, 4, avenue Blaise Pascal, TSA 60026, CS 60026, 63178 - Aubière Cedex, FRANCE <sup>7</sup> CERN, 1211 Geneva 23, Switzerland 8 Tomsk State University, 634050 Tomsk, Russia 9 Medical University of Vienna, Department of Radiation Oncology and Christian Doppler Laboratory for Medical Radiation Research for Radiation Oncology, Währinger Gürtel 18-20, 1090 Wien, Vienna, Austria <sup>10</sup> Université de Lyon, CREATIS; CNRS UMR5220, Inserm U1044, INSA-Lyon, Université Lyon 1, Centre Léon Bérard, France <sup>11</sup> Centre Antoine LACASSAGNE, Université Côte d'Azur – Fédération Claude Lalanne, Nice (France)

\*carlawinterhalter@gmail.com; adam.aitkenhead@christie.nhs.uk

# **Purpose**

This study investigates the influence of GATE/Geant4 settings on dose results and calculation times for clinical proton therapy.

A paper giving full details of the presented work is currently under review.

## **Materials and Methods**

- Same GATE/Geant4 versions as used in GATE-RTion1.0 [1]: GATE8.1 and Geant4.10.3.p03.
- Investigated physics lists (QGSP\_BIC\_EMY, QGSP\_BIC\_EMZ, QGSP\_BIC\_HP\_EMZ).
- Investigated production cuts on photons, electrons and positrons (restrict which particles are explicitly produced, 0.1mm-10mm).
- Simulations in a solid water phantom (for patient specific quality assurance) and within the patient CT (treatment sites: Neck, spine with titanium implant, paranasal), 0.25% statistical uncertainty at the 90-100% isodose level.

Figure 1: Calculation time & dose

difference due to physics settings

(phantom/range shifter, world).

and production cuts in


## **Conclusion**

#### **Physics lists:**

No substantial dose difference when using the HP (High Precision neutron) module. Best agreement to measurement data was found for QGSP BIC EMZ reference physics lists at the cost of 10-30% increased execution times compared to QGSP\_BIC\_EMY.

Production cuts on electrons, photons and positrons of 1 mm (phantom/rangeshifter) and 10 mm (world) substantially reduce calculation times.


# Comparison within solid water (patient specific quality assurance)



## **Physics lists:**

- · No substantial dose differences in the target volume when using the HP module, calculation times are increased by a factor of 1.1.
- Difference in scattering when comparing EMY against EMZ, calculation times are decreased by a factor of 1.3.

### **Production cuts:**



EMZ (0.1mm, 1mm) HP\_EMZ (0.1mm, 1mm) - EMY (0.1mm, 1mm) - EMZ (1mm, 10mm)

EMZ (0.1mm, 1mm)

# **Physics lists:**

- The high precision neutron libraries (HP) do not cause substantial dose differences in the field, calculation times are increased by a factor of 1.2.
- In agreement with the solid water calculations, there is a difference in scattering when comparing EMY against EMZ. Calculation times are decreased by a factor of 1.2.

#### **Production cuts:**

Dose differences at interfaces between high and low density tissues when comparing production cuts. Calculation times are decreased by a factor of 2.2.

Figure 2: Calculation time & dose difference in patient CT due to physics settings and production cuts in (phantom, world).

• No substantial influence on dose, calculation times are decreased by a factor of 4.7.







References [1] L. Grevillot et al., Technical Note: GATE-RTion: a GATE/Geant4 release for clinical applications in Scanned Ion Beam Therapy, accepted in Medical Physics (2020).

EMZ (0.1mm, 1mm) EMZ (0.1mm, 1mm)